

InteleConnect Quick Reference

Quick Tips

Logging in for the First Time via Computer






1. In your web browser, enter the web address for InteleConnect.

<https://pacs.premierradiology.com/Inteleconnect>

2. Do one of the following:
 - If you do not have an account for InteleConnect, click Apply for Account and follow the directions.
 - If your radiology provider has already given you an account, you can log in normally by typing your username and password, then pressing **ENTER**.

Logging in for the First Time via iPhone or iPad Device

1. From your app store, download the **InteleConnect app**.
2. Open the app and enter the following information:
 - Server: <https://pacs.premierradiology.com>
 - User Name: _____
 - Password: _____
3. Choose **My InteleConnect**
4. Log in with your user name and password.
5. These icons will allow you to do the following:

Click this:	To do this:
	Receive notifications
	Receive critical notifications
	Search by patient
	Change account preferences and notifications
	Log out

Setting Contact Information

1. Click Profile.
2. Click Account Preferences.
3. Type or edit your email address, telephone numbers, or address in the appropriate fields.
4. Click Save.

Setting Notification Preferences

1. Click Profile.
2. Click Notifications.
3. Enable or disable notifications.
4. Set the notification delivery option.
5. Click Save.

Viewing Notifications

1. In the Show Notifications list, choose one of the following options:

Choose this option:	To view:
All Notifications	All new and recent information sent to you from your radiology provider regarding your patients.
My Notifications	Only your preferred notifications; that is, the notifications you selected in your Profile.
Critical Notifications	Urgent or important information that requires your immediate attention and acknowledgement.
Pending Critical Notifications	Only pending critical results; that is, critical results that have not been completed or acknowledged.

2. Do one of the following:

To do this:	Click
See detailed information about the order	Status of the order for which you have been notified.
See detailed information about the order for which there is a critical result	Critical Result link.
See the list of the patient's current and prior orders	Name of the patient for which you have been notified.

3. To return to the list of notifications, click Notifications.

What's New!

- Patients by Location list (desktop version)
- My Patients list (desktop version)
- Window levelling and window level presets (desktop version)
- Access to orders that do not have studies









IntelConnect Quick Reference




Quick Tips

Viewing Images

Note: These images are for review only. They are not appropriate for diagnostic purposes.

- Do one of the following:
 - Search for a patient and then select the order you want to view.
 - Click the entry for the order in the Notification page.
 - Select a patient from the My Patient list or Patient by Location list and then select the order.
- Do one of the following:
 - To view thumbnails for all images for this order, click All Images.
 - To view thumbnails for key images for this order, click Key Images
- Click the thumbnail of the image you want to view.
- Do any of the following:

Click this:	To do this:
	Return to the Thumbnail view.
	View the previous series.
	View the next series.
	Scroll among the images in the series.
	Select the Window Levelling mode or the preset values. Click and drag the image up or down to adjust the brightness and left or right to adjust the contrast. Note: Window levelling is only available on desktop versions.
	Select the Pan mode. Click and drag the image in the required direction.
	Zoom in.
	Zoom out.

Click this:	To do this:
	Revert the image.
	Use the entire screen or browser window to view the image. You can also double-click the image.
	Exit the full screen. You can also double-click the image.

Printing Reports

- Click the entry for the order in the Notification page or search for the patient. The patient's report and thumbnails for key images appear.
- Click View PDF. The View PDF options dialog appears.
- Enable the checkboxes for the items that you want to appear in the printable report.
- Click View PDF.
- Use your browser's features to save or print the PDF.